Handi-Craft Company St. Louis, Missouri 63116 USA/EE.UU.
drbrownsbaby.com (800) 778-9001
Made in China. Fabricado en China.
May be patented in the USA, see drbrownsbaby.com/patents . Puede estar patentado en los EE.UU., vea drbrownsbaby.com/patents .
©2020 Handi-Craft Company. All Rights Reserved.
©2020 Handi-Craft Company. Todos los derechos reservados. F1

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Manufactured By:
Guangzhou Yongyi Industrial Co., Ltd.
B3 Factory, No. 29 Dongsan Street, Qixin Road,
Dongchong, Nansha Distric, GuangZhou, 511453
GuangDong, P.R.C.

DrBrown's

natural flow®

Single Electric Breast Pump

with SoftShape™ Silicone Shields

Model Numbers BF103-E/F-INTL and BF103-G-INTL

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Warning

- If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.
- Powered breast pumps that are designed for single users should never be rented or shared.
- This breast pump should not be operated while traveling on an airplane.

Contraindications

There are no known contraindications.

Precautions

- If discomfort is felt, suction can be broken by inserting a finger between the breast and the breast shield.
- If irritation or discomfort continues to occur, discontinue use and contact your healthcare provider.
- Before each use, visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.
- Inspect the tubing. If you find condensation in the air tube, wash and dry it immediately or replace the tubing. If you see breast milk inside the tubing or tubing connector, do not attempt to wash, clean or sterilize the tubing – contact Dr. Brown's Customer Service. Outside of the US and Canada, contact your local distributor.

Intended Use

The Dr. Brown's Natural Flow® Single Electric Breast Pump is a powered breast pump used by lactating women to express and collect milk from their breasts. This pump is intended for a single user.

Product Description

The Dr. Brown's Natural Flow® Single Electric Breast Pump is a personal use electric breast pump capable of single pumping. The device is electrically powered from an external AC powered supply and is provided non-sterile.

The Single Electric Breast Pump consists of the pump unit, breast shields, tubing, bottle vessel, duckbill valve, and other related accessories. The materials used are silicone (SoftShape™ Silicone Shields, nipples, membranes, duckbill valves), polypropylene (pump body and bottle vessels) and acrylonitrile butadiene styrene (pump unit). The electric breast pump's drive unit employs a diaphragm type vacuum pump, powered by a DC motor, supervised by a microcontroller. The microcontroller provides control over motor speed (vacuum creation) and solenoid (vacuum release).

The pump is reusable and intended for direct user contact in a home environment. The pump contains both backflow and overflow protection. The pump is not meant to be sterilized.

Additional parameters: Maximum suction pressure with no software control is -295mmHg. The shelf life of the pump is 2 years. This pump complies with safety standards IEC 60601-1 and IEC 60601-1-11.

Model BF103 has 3 modes:

1. Letdown Massage Mode – used to stimulate the release of milk from your breast. This mode has a quick cycle rate and helps prompt the reflex that happens when nerves in your breasts are stimulated (usually by your baby's sucking).
2. Expression Mode 1 – used after breast is stimulated to release milk at an intermediate rate.
3. Expression Mode 2 – used after breast is stimulated to release milk at a slow cycle rate.

PLEASE DO NOT RETURN THIS PUMP TO THE RETAILER

If you experience any issues with your Dr. Brown's Pump, please contact our dedicated Breast Pump Customer Service team by calling 1-833-224-7675.
Outside of the US and Canada, contact your local distributor.

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Table of Contents

Safety Information	3
Breast Pump Components.....	6
Cleaning.....	8
Assembly	10
Battery	14
Using the Pump	15
Pump Help	18
Storing & Feeding Breast Milk.....	19
Technical Data	21
Electromagnetic Emissions.....	22
Troubleshooting	24
Important Symbols.....	25
Disposal	26

Read all instructions before using this product.

WARNING: Read this ENTIRE Instruction Manual to become familiar with the features of the product before operating. Failure to operate the product correctly can result in damage to the product, personal property and cause serious injury. Failure to operate this product in a safe and responsible manner could result in injury or damage to the product or other property. This product is not intended for use by children. Do not use with incompatible components. This Manual contains instructions for safety, operation and maintenance. It is essential to read and follow all the instructions and warnings in the manual, prior to setup or use, in order to operate correctly and avoid damage or serious injury.

DANGER - To reduce the risk of electrocution:

- Always unplug electrical product immediately after use.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place or drop into water or other liquid.
- Do not reach for a product if it has fallen into water. Immediately unplug from electrical outlet.

WARNING - To avoid fire, electrocution, or injury:

- Danger of electrocution! Keep the device dry! Never immerse in water or other liquids! Do not run water over the pump.
- The breast pump and accessories are not heat-resistant: keep away from radiators, open flame and sunlight.
- Keep the power adapter away from heated surfaces.
- Do not reach for any electrical device if it has fallen into water. Unplug immediately.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- The breast pump should never be left unattended when plugged into a power source.

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Read all instructions before using this product.

WARNING (Cont):

- Never operate an electrical device if it has a damaged cord, plug or battery pack, if it is not working properly, if it has been dropped or damaged, or dropped into water. If damage is found, immediately discontinue use of power adapter or battery pack and contact your local distributor.
- Do not use the breast pump while bathing or showering.
- Close supervision is necessary when the breast pump is used in the vicinity of children to prevent strangulation due to power adapter cord or tubing.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and under the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- Follow the exact order of the steps in powering breast pump as listed on page 13 of this Instruction Manual.
- For the purpose of recharging the battery only use the power adapter supplied with the breast pump.
- Make sure the voltage of the power adapter is compatible with the power source.
- This appliance contains batteries that are non-replaceable.

CAUTION - To avoid risk of minor injury:

- Wash hands thoroughly with soap and water before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.
- Separate and wash all parts that come in contact with breast milk immediately after use to avoid dry up of breast milk residue and prevent growth of bacteria. Always inspect SoftShape™ Silicone Shield, connector, valve, membrane and tubing prior to use for cleanliness.
- Only use drinking-quality tap or bottled water for cleaning your breast pump and parts.
- Use the breast pump only for its intended use as described in this manual.
- Do not use the breast pump while sleeping or overly drowsy.
- This is a single user product. Use by more than one person may present a health risk.

Read all instructions before using this product.

CAUTION (Cont):

- Do not drive while pumping.
- Never use while pregnant, as pumping can induce labor.
- Clean and sanitize all parts that come in contact with your breast and breast milk prior to first use.
- Wash all parts that come in contact with your breast and breast milk after every use.
- Do not store wet or damp parts, as mold may develop.
- Inspect all appropriate pump components before each use.
- Do not run pump with wet tubing. Doing so will harm the pump motor.
- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain.
- If you are experiencing discomfort at the base of the nipple due to rubbing of your breast tissue with the breast shield tunnel, discontinue use and contact your healthcare professional or lactation consultant.
- For assistance with correct breast shield sizing and comfort, please contact your local distributor or contact a lactation consultant.
- Contact your healthcare professional or lactation consultant if you can express only minimal or no milk or if expression is painful.
- Do not try to express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.
- Do not thaw frozen breast milk in a microwave or in a pan of boiling water.
- Do not microwave breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breast milk.)

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
Cleaning

BEFORE FIRST USE, AND AFTER EVERY USE

Clean and sterilize all pump parts that come in contact with milk:

ALWAYS USE DRINKING QUALITY WATER FOR CLEANING!

Cleaning Instructions for Milk-Contact Parts:


1. Disassemble pump following instructions on page 10.
2. Rinse milk residue from parts under cool running water.
3. Do not place parts directly in the sink. Soak in warm soapy water for several minutes.
4. Use a cloth or soft brush used only for cleaning infant feeding items to clean all parts.
5. Rinse parts under hot running water until all soap has been removed, being careful that temperature of water will not burn or irritate skin.
6. Alternately, wash all parts in the top rack of the dishwasher. Place small items in a closed-top basket or mesh bag. Run dishwasher using hot water and a heated drying/sanitizing cycle.
7. Rest parts on a clean paper towel or rack to air dry. All parts must be completely dry before storing to prevent bacteria and mild growth. Store clean, dry parts in a bag or container until next use.
8. Sanitize/sterilize before each use. (Microwave and electric steam sterilizers may be used to sterilize milk contact parts.)
9. For convenience, special wipes for cleaning breast pump parts may be used when away from home. However, all parts must be cleaned using soapy water or in the dishwasher before pumping again.

NOTE

- Please clean all accessories according to cleaning guide every time before using.
- Let parts cool after washing in dishwasher or sterilizing. Handle carefully.
- Avoid using cloth towels to dry your pump parts because they can carry germs and bacteria that are harmful to your breast milk and your baby.
- Try not to touch the inside of any parts that will come in contact with your breast milk.
- See separate instructions for bottle cleaning & assembly.

NON MILK-CONTACT PARTS:

The following parts do not come in contact milk and should only be cleaned as needed.

Special Cleaning Instructions for Tubing:

1. Tubing does not come in contact with milk and should only be washed if dirty or condensation is present.
2. If condensation is from a previous washing process, there is no need to wash again.
3. Wash with warm soapy water and rinse clean with clear water.
4. Allow water to drain from tubing completely after washing.
5. Shake out any remaining water droplets and hang to dry.
6. Do not use tubes until completely dry.*

*A small amount of condensation is allowable; however, a tube with water droplets may damage the pump motor.

**DO NOT SUBMERGE THE
PUMP MOTOR IN WATER!
WIPE CLEAN WITH
DAMP CLOTH ONLY
AFTER EACH USE.**

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

DISASSEMBLY:

The SoftShape™ Silicone Shield/Milk Collection Kit comes pre-assembled. All milk-contact components must be disassembled, washed and sterilized prior to first use.


Wash hands before touching & assembling milk collection kit.

All milk-contact components must be disassembled, washed and sterilized prior to first use.

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


TUBING:



POWER ADAPTER:

To prevent damage to your power adapter cord, always grasp the body of the adapter as shown above.
Never pull by the cord to unplug.
Always unplug the adapter from the pump before storage.

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
Battery

BATTERY SPECIFICATION:


Dr. Brown's Natural Flow® Single Electric Pump has a built-in rechargeable battery for your convenience. A single charge of the battery lasts approximately 1.5 hours, or about 6 pumping sessions of 15 minutes each.

BATTERY LIFE INDICATOR:

When using your pump while it is not plugged in, check the Battery Level Indicator on the display screen to see how full the battery is.

To charge your pump's battery, simply plug in the included power adapter as shown on page 13. While charging, the display will show both the Power Indicator and Battery Level Indicator. The battery is fully charged when the Battery Level Indicator shows four solid bars.


Using the Pump

Before pumping, wash hands thoroughly with warm soapy water. Inspect all breast pump parts for cleanliness - IF UNWASHED, FOLLOW INSTRUCTIONS ON PAGE 8. Look closely at the silicone shield, membrane, and duckbill valve to inspect for any signs of wear or damage. Suction may decrease if parts are damaged or not properly assembled. Contact your local distributor to inquire about replacement parts.

Do not try to express milk if the level of vacuum is painful, this may lead to nipple and/or breast trauma and cause difficulties with expression. Contact your medical professional or lactation consultant if you can only express minimal amounts of milk, or if expression is painful in any way.

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




LETDOWN and EXPRESSION MODE:

Letdown Mode is used to stimulate the release of milk from the breast. This mode mimics baby's way of rapidly sucking at the beginning of a breastfeeding session, prompting the reflex that releases breast milk.

Expression Modes 1 and 2 are used once the breast milk has been released and is now expressing from the nipple. This mode mimics baby's more relaxed feeding that happens once milk is flowing and draws milk more deeply from the breast. Expression Mode 1 operates at an intermediate cycle rate, and Expression Mode 2 operates at a slow cycle rate.

Using Letdown Mode:

1. Position the assembled milk collection kit on breast as directed on page 17.
2. Press the POWER BUTTON to turn the unit on.
3. The pump will automatically display Letdown Massage Mode at the last vacuum level used and the timer will begin.
4. If needed, adjust vacuum level to your preferred starting level using the VACUUM UP and VACUUM DOWN BUTTONS. There are five user-selectable vacuum levels.
5. Letdown Massage Mode will automatically change to Expression Mode 1 after two minutes. See Mode Elapse Time on display. To switch to Expression Mode 1 sooner, press the MODE SELECTOR BUTTON.

It may take several pumping sessions to discover your most comfortable and productive vacuum level.

Once letdown has begun (milk is visibly expressing from nipples), switch to Expression Mode.

Using Expression Mode 1:

1. Press the MODE SELECTOR BUTTON. The Expression Mode 1 Indicator Symbol will appear. The pump will automatically operate Expression Mode 1 at the last vacuum level used.
2. If needed, adjust vacuum level to your preferred level using the VACUUM UP and VACUUM DOWN BUTTONS. There are eight user-selectable vacuum levels.
3. To return to Letdown Massage Mode, press the MODE SELECTOR BUTTON.


Using Expression Mode 2:

1. To switch to Expression Mode 2 from any mode, press the MODE SELECTOR BUTTON and hold for two seconds. The Expression Mode 2 Indicator Symbol will appear. The pump will automatically operate Expression Mode 2 at the last vacuum level used.
2. If needed, adjust vacuum level to your preferred level using the VACUUM UP and VACUUM DOWN BUTTONS. There are eight user-selectable vacuum levels.

It may take several pumping sessions to discover your most comfortable and productive cycle and vacuum levels.

Placement of the Milk Collection Kit:

Position the SoftShape™ Silicone Shield over your breast, centering the nipple in the neck of the shield. Your nipple should not touch or rub on the side walls of the shield. Check nipple positioning again once the pump is operating, to ensure that it remains away from the sides. If not, try repositioning. Your breast should completely fill the SoftShape™ Silicone Shield so that a seal may form to achieve proper vacuum levels.

While you may feel some discomfort when first pumping, it should not cause pain. Adjust the vacuum levels to lower settings to achieve your preferred comfort level and milk flow.

TIPS FOR EXPRESSING BREAST MILK

Relax! Create a comfortable and relaxing setting for your pumping sessions.


Stay Hydrated. Drink plenty of water throughout the day and keep a glass of water with you while you pump.

Create a Routine. Condition yourself to express milk at the same time each day by creating a routine around pumping. Try pumping each morning after you shower or pumping in the afternoon while watching your favorite television show.

Think of Baby! Simply imagining that you are breastfeeding may help induce letdown. If you are away from home, some moms find that a picture of baby or a recording of their cry helps.


Massage. A tried and true method - massaging your breasts before you begin pumping can help to quicken letdown. Massage again when you see your milk slowing to help release all milk from the breast.

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COMMONLY ASKED QUESTIONS

What is the difference between Letdown Mode and Expression Mode?

Letdown Mode helps stimulate letdown (milk releasing from breast). Expression Mode is used after Letdown Mode - once the milk is visibly flowing from nipple - to draw milk from breast for collection.

What is the difference between cycle and vacuum levels?

Vacuum levels represent the amount of suction. Cycle levels represent the speed of the on/off suction cycle.

How do I find my preferred cycle and vacuum levels?

Every mom reacts differently to cycle and vacuum levels. Start low and work your way up to a comfortable and productive level. It may take a few pumping sessions to find your ideal settings, and your ideal settings may change over time. Dr. Brown's™ Single Electric Pump settings are customizable so that you can have your perfect pumping experience every time!

What do I do if milk gets in the tubing?

Wash the tubes as directed on page 9 of this instruction manual. Remember to dry the tubing completely before using - water in the tubes may damage the motor. If mold appears in tubing, discard and replace.

What do I do if there is low/no suction?

Check that the tubing is securely connected. Disassemble and check parts for damage. Reassemble, ensuring that parts are connected properly, creating a good seal for the vacuum. While pumping, ensure that your breast is also creating a seal by fully filling the shield.

Why does my pump turn off after I have been pumping for a while?

Your pump has a safety feature that shuts off the pump after 30 minutes.

If I have a problem, can I take my pump to a local repair center?

No. Please contact your local distributor for more information.

After every use

- Store milk safely. Seal milk collection bottle or milk collection bag, label with date and time, and immediately place in a refrigerator, freezer, or cooler bag with ice packs. Refer to "Human Milk Storage Guidelines" on page 19 for more details.
- Wipe down breast pump with a clean soft towel, and disinfect pumping area.
- Take apart breast pump tubing and separate all parts that come in contact with breast/breast milk.
- Rinse breast pump parts that come into contact with breast/breast milk by holding under running water to remove remaining milk.
- Clean pump parts that come into contact with breast/breast milk as soon as possible after pumping. Follow the cleaning steps provided on page 8.

Guidelines Are for Home Use and Not for Hospital Use.

Use the following guidelines for the storage of breast milk.

By following recommended storage and preparation techniques, breastfeeding mothers and caretakers of breastfed infants and children can maintain the safety and quality of expressed breast milk for the health of the baby. These are general guidelines for storing breast milk at different temperatures. Various factors (milk volume, room temperature when milk is expressed, temperature fluctuations in the refrigerator and freezer, and cleanliness of the environment) can affect how long breast milk can be stored safely. Always keep milk covered when storing.


Human Milk Storage Guidelines			
Storage Locations and Temperatures			
Type of Breast Milk	Countertop 77°F (25°C) or colder (room temperature)	Refrigerator 40°F (4°C)	Freezer 0°F (-18°C)
Freshly Expressed or Pumped	Up to 4 hours	Up to 4 days	Within 6 months is best Up to 12 months is acceptable
Thawed or Previously Pumped	1-2 hours	Up to 1 day (24 hours)	Never refreeze human milk after it has been thawed
Leftover from a Feeding (baby did not finish bottle)	Use within 2 hours after the baby is finished feeding		

Adapted from "ABM Clinical Protocol #8: Human Milk Storage Information for Home Use For Full-Term Infants," Revised 2017

Breast Milk Storage Guidelines Center for Disease Control and Prevention ("CDC");
https://www.cdc.gov/breastfeeding/recommendations/handling_breastmilk.htm

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
STORE

- Store breast milk in storage bags made for freezing human milk or clean glass or hard BPA-free plastic bottles with tight fitting lids. Do not use containers with the number 7 which may contain BPA. Do not use disposable bottle liners or other plastic bags to store breast milk.
- Label milk with date and baby's name if delivering to childcare.
- After each pumping session you can:
 - Keep milk at room temperature. Breast milk is safe for up to 4 hours at room temperature (up to 77°F or 25° C) after pumping.
 - Refrigerate it. Breast milk is safe in the refrigerator for up to four days.
 - Place milk in the freezer. If you do not plan to use milk within four days of pumping, freeze right away.
 - Use cooler packs. You can store breast milk in an insulated cooler with frozen ice packs for up to 24 hours after pumping. Store in refrigerator or freezer after 24 hours.
 - Do not store milk on the door of the refrigerator or freezer.
- Freeze milk in small amounts of 2 to 4 ounces (59 to 118 milliliters) to avoid waste when warming.
- Breastmilk expands when frozen. Leave a little space at the top of your container. Wait to tighten caps until milk is completely frozen.

THAW

- Always thaw and use the oldest milk first.
- Thaw milk under warm running water, in a container of warm water, in the refrigerator overnight, or in breast milk warming machine.
- Never thaw milk in the microwave. Microwaving destroys nutrients baby needs, and hot spots created in the microwave can burn baby's mouth.
- Use milk within 24 hours of thawing.
- Use thawed or warmed milk within two hours.
- Never refreeze thawed milk.


FEED

- To heat milk, place the sealed container under warm running water, in a container of warm water, or in breast milk warming machine
- Do not heat milk on the stove or in the microwave.
- Test the temperature before feeding it to your baby by putting a few drops on your wrist.
- Swirl the milk, mixing the fat back into the milk.
- Leftover milk should not be used after two hours.


Single Electric Pump Motor

Vacuum Range: -35 to -230 mmHg
Cycle Rate: 40 – 63 cpm
Mode of operation: Continuous
Weight: 340 g

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

Electromagnetic Emissions

The Dr. Brown's Natural Flow® Single Electric Breast Pump is intended for use in the electromagnetic environment specified below. The user of the Breast Pump should assure that it is used in such an environment.

Guidance and manufacturer's declaration – electromagnetic emissions		
Emissions	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Breast Pump uses pull of vacuum only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Breast Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	± 1 kV line(s) and neutral	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5s	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Breast Pump requires continued operation during power mains interruptions, it is recommended that the Breast Pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the AC mains voltage prior to application of the test level


Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3V rms 150 kHz to 80MHz	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the Breast Pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 1,2\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation Distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of Equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Breast Pump is used exceeds the applicable RF compliance level above, the Electric Breast Pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Breast Pump.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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ISSUE	POSSIBLE CAUSES	CHECKS	SOLUTION
No Power	Adapter	Check if the power cord is not connected or it is damaged.	Confirm if plugged into working wall outlet and in power port in back of the pump.
	Battery Dead	Plug power adapter into wall and pump motor. Screen should light up and battery level indicator should appear.	Allow battery to charge or use pump while it is plugged into an outlet.
	Battery Malfunction	Perform checks for "Battery Dead." If these do not fix issue, it is possible the battery has reached its end of life (approximately 2 years).	Consider purchasing new device, or if pump is still under warranty, contact local distributor.
No Suction	Tubing	Check connection to membrane cover and to pump port on the back of the pump motor.	Confirm connections.
Low Suction	Tubing	Inspect tubing for holes or cracks.	Replace tubing if crack is in the middle of the length. If cracked near tube ends, remove end and trim away damaged tube and reinstall tube end.
	Assembly error	Check for loose assembly of milk collection kit components.	Disassemble and reassemble. Take care to tighten the membrane cap.
	Duckbill valve	Check for damage.	Replace valve if damaged.
	SoftShape™ Silicone Shield	Check for damage. Softshape™ Silicone Shield is completely on the breast.	Replace if damaged. Press more firmly against the breast.
	Silicone membrane	Check for damage.	Replace membrane if damaged.
Pain	Too high suction	Check the vacuum level.	Set to the lowest suction level or consult your lactation consultant or medical professional.

	This symbol indicates do not dispose the device together with unsorted municipal waste.
	This symbol indicates compliance with international requirements for protection from electric shock.
	This symbol indicates that the power adapter is a class II device.
	This symbol indicates that the power adapter is for indoor use only.
	Certification mark that indicates conformity with health, safety, and environmental protection standards for products sold within the European Economic Area.
	This symbol indicates polarity of D.C. power connector.
	This symbol indicates alternating current.
	This symbol indicates direct current.
	This symbol indicates that this device should be kept dry.
	This symbol indicates the need for the user to consult the instructions for use.
IP22	This symbol indicates the device meets IP22 classification.
	This symbol indicates when the device was manufactured.
	This symbol indicates where the device was manufactured.
SN	This symbol indicates the manufacturer's serial number that identifies the device.
LOT	This symbol indicates the manufacturer's batch code so that the batch or lot can be identified.
REF	This symbol indicates the manufacturer's catalogue number.
	This symbol indicates that a specific separate supply unit is required for connecting electrical equipment to the supply mains.

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Disposal of Breast Pump

The breast pump unit is made of a combination of plastic and metal, and contains a lithium battery. Dispose in the manner required by your local, state, and national laws.

Disposal of Collection Kit Items

The hard, plastic pieces of the collection kit can be recycled per your local regulations. This includes the following parts, as outlined on page 6.

- Membrane Cap
- Pump Connector
- Bottle Vessel

The soft, silicone pieces of the collection kit cannot be recycled. This includes the following parts, as outlined on page 6.

- SoftShape™ Silicone Shield
- Membrane
- Tubing
- Duckbill Valve

Disposal of Dr. Brown's® Options+™ Bottle and Accessories

The hard, plastic accessories and pieces of the bottle can be recycled per your local regulations. This includes the following parts, as outlined on page 6 and in the included separate bottle instructions.

- Caps
- Collars
- Reservoir
- Storage/Travel Cap
- Narrow-to-Wide Neck Adapter

The soft, silicone pieces of the bottle and the bottle cleaning brush cannot be recycled. This includes the following parts, as outlined on page 6 and in the included separate bottle instructions.

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DR. BROWN'S and all other trademarks are owned or used under license by Handi-Craft Company.